## INSTRUCTIONS FOR USE

### INTENDED USE

ESR-Chex® Plus is an assayed bi-level control for evaluating the accuracy and precision of automated Erythrocyte Sedimentation Rate (ESR) methods.

### SUMMARY

Clinical laboratories are required to establish quality control programs for automated and semi-automated procedures. Daily use of ESR-Chex Plus provides assurance that the erythrocyte sedimentation rate procedure is being performed properly.

### PRINCIPLES

ESR-Chex Plus is designed to be handled in the same manner as a patient specimen. There are several factors concerning the environment and techniques that may affect the outcome of erythrocyte sedimentation rate methods. When these are properly controlled, the values recovered with ESR-Chex Plus will be within the Expected Ranges printed on the assay. Accurate analysis of patient specimens is assured.

### REAGENTS

ESR-Chex Plus contains stabilized human red blood cells in a preservative medium.

### PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
5. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

ESR-Chex Plus is stable through the expiration date when stored at 2 °C to 10 °C. After bringing to room temperature, ESR-Chex Plus is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 18 °C to 30 °C. Open-vial stability for this control is defined as the maximum number of continuous days a vial can be brought to room temperature and mixed for analysis. Vials must be discarded at the lot expiration or the OVS expiration, whichever comes first. Working vials are to be stored at room temperature, and may not be returned to refrigerated temperature for additional use.

### INDICATIONS OF PRODUCT DETERIORATION

Discoloration of the product may be caused by overheating or freezing during shipping or storage. Gross hemolysis (darkly colored supernatant) may be indicative of product deterioration. However, moderately colored supernatant is normal and should not be confused with deterioration of the product.

### INSTRUCTIONS FOR USE

Follow the instrument manufacturer's instructions for ESR testing for automated methods.

#### MINI-CUBE Instrument Procedure

Utilize the MINI-CUBE Barcodes on the assay sheet to register each new QC lot and perform QC analysis. DO NOT use the barcodes on the control vials, as the vial barcodes are exclusively for use with the CUBE 30 Touch.

Note: Failure to scan the correct barcode may result in the incorrect assay range applied to the result. Insert the control tube with the label-free gap facing the dot on the right side of each well.

#### CUBE 30 Touch Instrument Procedure

ESR-Chex Plus controls are automatically registered and processed in a QC file when the barcode is scanned. Position the sample so the barcode faces the mark (|) found at each tube position on the carousel.

### MIXING AND HANDLING INSTRUCTIONS

1. Remove vials from refrigerator and allow them to equilibrate to room temperature (15-20 minutes).
2. To Mix:

For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).

- a. Vortex unopened vials for up to 60 seconds to resuspend product prior to first time use.
- b. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- c. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- d. Examine the bottom of the vial to ensure that the product is completely resuspended and repeat steps 2b-2c, if necessary.
  - e. Analyze the sample immediately after mixing. If mixed vials sit for more than 1 minute before analysis, remix the vial by inverting 8-10 times.
3. Working vials stored at room temperature (18 °C to 30 °C) should not require additional vortexing. Mix working vials by following steps 2b-2e and repeat as necessary to ensure the cells are completely resuspended.
  4. It is recommended that working vials be stored at room temperature (18 °C to 30 °C) for the duration of the open-vial stability and subsequent analysis.
  5. Follow the manufacturer's directions for the method in use to perform sedimentation rate testing.

### LIMITATIONS

ESR-Chex Plus is to be used for Erythrocyte Sedimentation Rate testing only. It is not intended for controlling other hematology procedures.

### EXPECTED RESULTS

The assay values are derived from replicate analysis on automated methods. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits. However, the control means established by the laboratory should fall within the Expected Range specified for the control.

In order to achieve values within the published assay ranges, laboratories should not use a working vial of ESR-Chex Plus longer than the open-vial dating, as indicated on the assay sheet. Failure to adhere to this guideline may result in values outside the published expected ranges. It is advised that laboratories establish their own means and ranges when utilizing control material outside of these handling instructions.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2<sup>1</sup> it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

If established means do not fall within the Expected Range, please evaluate all mechanical and physical factors that could affect the outcome, such as handling, temperature, vibration, tube position and product expiration. Inadequate mixing of the product will affect the outcome of the test, see Handling Instructions. Sample must be analyzed within 60 seconds of thoroughly mixing the vial. Assay a new working vial of ESR-Chex Plus. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

### REFERENCES

1. Clinical and Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.

### QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or [statsdata@streck.com](mailto:statsdata@streck.com). Additional information can be found at [streck.com](http://streck.com).

### ORDERING INFORMATION

Please call our Customer Service Department 800-228-6090 for assistance. Additional information can be found online at [streck.com](http://streck.com).

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com).

See [streck.com/patents](http://streck.com/patents) for patents that may be applicable to this product.

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#### Procedure for MINI-CUBE instrument

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#### Procedure for CUBE 30 Touch instrument

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